

# Laboratory Muffle Furnace

Laboratory muffle furnace for ashing, analysis and small-part heat treatment. Match the process to the 1000 °C rated maximum and chamber dimensions.



## Technical Information

|                            |                    |
|----------------------------|--------------------|
| Rated temperature          | 1000 °C            |
| Rated power                | 2.5 kW             |
| Rated voltage              | 220 V              |
| Working voltage            | 220 V              |
| Empty-chamber heating time | ≤60 min            |
| Chamber dimensions         | 200 × 120 × 80 mm  |
| Dimensions                 | 595 × 410 × 467 mm |
| Gross weight               | 80 kg              |

## Standard Configuration

The quotation must identify the selected instrument configuration and supplied accessories. Where several specifications are listed, they describe separate versions in the same order; they do not describe one combined configuration.

## Optional Configuration

Optional equipment and consumables are supplied only when itemized in the quotation. Confirm dimensions, compatibility, electrical supply and the required configuration before ordering.

---

## Selection and ordering

Include the PRECETTA model, selected specification, sample or application, quantity and destination in your inquiry. Use the complete PDF to compare variants and configuration details.

## Technical Information

Specifications apply to the stated chamber and rating. Confirm the controller and power connection. Crucibles and sample holders are not stated as included.