

Sulfur Dioxide Distillation System Series

Multi-position distillation for sulfur dioxide residue sample preparation in food and herbal materials. Select the position count and flask-capacity combination.



Technical Information

Distillation positions	4 / 6 / 6
Distillation flask capacity	1000 / 500 / 1000 mL
Receiving flasks	4 × 100 / 6 × 100 / 6 × 100 mL
Electrical supply	220 V; 50 Hz
Temperature limit	up to 200 °C
Display	touchscreen
Dry-boil protection	yes
Leakage protection	yes
Cooling circulation	integrated pump

Standard Configuration

The quotation must identify the selected instrument configuration and supplied accessories. Where several specifications are listed, they describe separate versions in the same order; they do not describe one combined configuration.

Optional Configuration

Optional equipment and consumables are supplied only when itemized in the quotation. Confirm dimensions, compatibility, electrical supply and the required configuration before ordering.

Selection and ordering

Include the PRECETTA model, selected specification, sample or application, quantity and destination in your inquiry. Use the complete PDF to compare variants and configuration details.

Technical Information

The three versions are 4 × 1000 mL, 6 × 500 mL and 6 × 1000 mL distillation arrangements. Nitrogen generators or cylinders and external cooling units require separate selection. Confirm the integrated circulation pump and any external cooling equipment in the delivery list.