
Endotoxin Test Tube Adapter

Endotoxin Test Tube Adapter. Replacement accessory for the laboratory instrument configuration specified below. Use the selected instrument model and the existing component dimensions or connector details to confirm the correct replacement.



Key features

- Endotoxin Test Tube Adapter
- T-END-002

Technical Information

Accessory model	T-AEN-005
Minimum order quantity	1
Documented host model	T-END-002

Standard Configuration

Supply of the selected accessory, in the confirmed order quantity. The host instrument and other accessories are excluded unless itemized in the quotation.

Optional Accessories

Related replacement parts are quoted separately after compatibility review. No unlisted fitting or installation service is assumed included.

Selection / Ordering Information

Include the PRECETTA instrument model, serial number, component photograph and required quantity in the inquiry. The photograph identifies the accessory type; confirm the exact ordered variant against its specification before procurement.