

Aqueous Sample Syringe Filters

Aqueous Sample Syringe Filters. Replacement accessory for the documented laboratory application. Compare the specification options below and confirm the instrument configuration before selecting a variant.



Key features

- 0.22 / 0.45 / 0.8 μm
- Minimum order quantity: 100

Technical Information

Accessory model	T-ASM-002
Specification options	0.22 / 0.45 / 0.8 μm
Minimum order quantity	100

Standard Configuration

Supply of the selected accessory, in the confirmed order quantity. The host instrument and other accessories are excluded unless itemized in the quotation.

Optional Accessories

Related replacement parts are quoted separately after compatibility review. No unlisted fitting or installation service is assumed included.

Selection / Ordering Information

Include the PRECETTA instrument model, serial number, component photograph and required quantity in the inquiry. The photograph identifies the accessory type; confirm the exact ordered variant against its specification before procurement.