

Advanced Pyrogen Temperature Monitor

Records probe temperatures for grouped rabbit pyrogen test workflows. Configurable probe channels, test records and documented audit functions support laboratory review; confirm the supplied software and protocol.



Key features

- Interchangeable temperature probes
- Grouped test records and PDF reports
- Documented access controls and audit logs
- Battery operation depends on probe load

Technical Information

Temperature range	0-50 °C
Temperature accuracy at 32-44 °C	$\leq \pm 0.1$ °C
Probe capacity / standard supply	≤ 96 / 33
Probes per group	3
Touch display	12.1 in
Battery operation: probe dependent	5-10 h

Standard Configuration

Monitor with 33 probes stated as standard; expandable to 96 probes. Confirm room-monitoring, data-entry and wireless alarm modules and quantities in the quotation.

Optional Accessories

Calibration materials, replacement consumables and qualification documents are separately confirmed and quoted where available.

Selection / Ordering Information

Confirm probe count, grouping protocol, calibration intervals, required records and software version. Battery runtime and wireless alarm range depend on operating conditions. Audit functions require laboratory validation; they do not establish certified regulatory compliance.