

High-Speed Laboratory Mill

High-speed mill for dry laboratory materials in a sealed stainless steel cup. Match the sample mass and target fineness to the application.



Technical Information

Capacity	200 g
Fineness	30-200 mesh
Motor power	1000 W
Motor speed	25000 rpm
Grinding cup dimensions	135 × 100 mm
Dimensions	200 × 220 × 280 mm
Cup material	stainless steel
Typical grinding time	0.5-3 min

Standard Configuration

The quotation must identify the selected instrument configuration and supplied accessories. Where several specifications are listed, they describe separate versions in the same order; they do not describe one combined configuration.

Optional Configuration

Optional equipment and consumables are supplied only when itemized in the quotation. Confirm dimensions, compatibility, electrical supply and the required configuration before ordering.

Selection and ordering

Include the PRECETTA model, selected specification, sample or application, quantity and destination in your inquiry. Use the complete PDF to compare variants and configuration details.

Technical Information

The 200 g grinding cup is stainless steel, with a sealed quick-opening cover. The stated fineness and grinding time depend on the dry sample; do not combine this product with unspecified larger-capacity models.